

# SmartVue® Shielded Window System

**ADDING NATURAL LIGHT TO TREATMENT ROOMS**



## **veritas** **SmartVue®**

**BUILDING A BETTER PATIENT EXPERIENCE**

**Introducing SmartVue, an innovative new radiation treatment experience that eliminates the claustrophobic confines of the typical treatment room by bringing the great outdoors into your linac room.**

SmartVue is a breathtaking new advancement in patient treatment that allows natural light directly into the treatment room, providing the patient with a beautiful, uplifting view of the outside before and after the treatment process. This groundbreaking new shielded scenic wall system will add an outside view to your treatment room without affecting the integrity of the surrounding radiation shielding.





## HOW IT WORKS

As patients enter the room, they are greeted with a beautiful view to the outside as the room is bathed with natural light. As the treatment begins, the lights are dimmed and the glass of the window “frosts over,” hiding radiation shielding that moves in behind the glass. When the treatment is finished and the lights in the room come back up, the shielding can then slide out of sight as the glass clears to once again reveal the beautiful scenery outside the window.



### A Better Experience

Natural light positively affects how patients react and how their bodies respond to treatment. Bringing an outside view into the treatment room will not only positively influence the patient, but the healthcare personnel as well. You owe it to all of them to provide the most positive radiation treatment experience there is.

### Scenic Packages for Varying Site Conditions

Fully adaptable to site conditions, the SmartVue window wall is available in a variety of sizes and configurations. Optional enhancements to the patient experience are available through unique audio/lighting “experience” packages, variable window opacity and customizable graphics.

Veritas continues to lead the way to a better radiation shielding experience keeping pace with the remarkable advancements of the treatments themselves.



For more information on Veritas SmartVue Shielded Window Systems or how Veritas Medical Solutions can benefit your next project, contact us at **888.242.6760** or visit **[veritas-medicalsolutions.com](http://veritas-medicalsolutions.com)**